

Wavelength® Product Comparison Chart

	Wavelength LV (Low Viscosity) Ultrasound Gel	Wavelength MP (Multi-Purpose) Ultrasound Gel	Wavelength CL (Clear) Ultrasound Gel	Wavelength TMS Gel	Wavelength HV (High Viscosity) Gel	Wavelength Ultrasound Lotion
Colour	Green	Blue	Colourless	Slight yellow tint	Colourless	White
Viscosity @ 25°C	Low (40–50 kcps)	Medium (60–80 kcps)	Medium (60–80 kcps)	Very High (150–160 kcps)	High (80–110 kcps)	Low (20-40 kcps)
pH @ 25°C	Neutral 6.5 – 7.5	Neutral 6.5 – 7.5	Neutral 6.5 – 7.5	Neutral 6.5 – 7.5	Neutral 6.5 – 7.5	Neutral 6.5 – 7.5
Stability	Bacteriostatic	Bacteriostatic	Bacteriostatic	Bacteriostatic	Bacteriostatic	Bacteriostatic
Primary Application	Imaging Ultrasound	Therapeutic Ultrasound	Therapeutic Ultrasound	TENS Applications	Therapeutic Ultrasound	Therapeutic Ultrasound